



# REEVUE

## REEVUE Indirect Calorimetry

The REEVUE measures the oxygen that the body consumes. Using this measurement it calculates a patient's Resting Energy Expenditure (REE). This technique is referred to as Indirect Calorimetry and is used for nutritional assessment and medical nutrition therapy.

The REEVUE performs a simple 10 minute breath test. It displays results and provides a printout that is ideal for maintaining records, as well as helping the clinician and patient to understand the results.

The REEVUE is lightweight and portable. It does not require a computer or manual calibration.

The REEVUE provides an accurate method for measuring VO<sub>2</sub> (oxygen consumption) and REE (resting energy expenditure).

### Indications For Use

The REEVUE can be used for any patient where a measurement of metabolic rate or oxygen uptake (VO<sub>2</sub>) will be useful, including:

- weight management
- nutritional assessment services

**Call OPS Medical Today!**

**1 (800) 430-0048 or**

**email: [sales@opsmedical.com](mailto:sales@opsmedical.com)**



## KORR

www.opsmedical.com 1 (800) 430-0048



### Features

#### ACCURATE

The REEVUE has been validated against the clinical standard metabolic cart and has demonstrated equivalent accuracy in resting volunteers for both oxygen consumption and resting energy expenditure.

#### EASY TO USE

When the REEVUE is turned on, it automatically self-calibrates and is ready to use. After the test is completed, results are presented on the display screen and are printed directly to an optional printer. No computer or software is needed to operate the REEVUE.


#### AFFORDABLE

The REEVUE is priced at only a fraction of the cost of a standard metabolic cart.

# Clinical Validation of the REEVUE

A study comparing resting energy expenditure (REE) and oxygen consumption (VO<sub>2</sub>) measurements reported by the **REEVUE** and by the Deltatrac metabolic cart was done at the University of Utah. This study, which has been submitted for publication, shows very good agreement between the **REEVUE** and the Deltatrac for measurement of both VO<sub>2</sub> and REE.

**REEVUE vs. Deltatrac XY Plot, Energy Expenditure**



### Key Statistics:

Average difference:	24.7 kcal/day (-1.45%)
Standard deviation of difference:	57.1 kcal/day (2.85% of reading)
Correlation:	$r^2 = 0.975$ ( $r = 0.987$ )

### REEVUE Specifications

<b>Air Flow Sensor</b>	Resolution Accuracy Flow Range Breathing Rate Tidal Volume	10 ml / sec ±2% of reading -40 to 150 LPM 5 to 40 breaths/minute 200 to 3000 ml	Weight Measurement Time	5.75 lbs. (2.6kg) 10 minutes
<b>Oxygen Sensor</b>	Accuracy Resolution	±0.2 %O <sub>2</sub> 0.01 %O <sub>2</sub>	<b>Operating Environment</b>	Temperature Range Barometric Pressure Range Relative Humidity Range
<b>Measurement Module</b>	VO <sub>2</sub> Range VO <sub>2</sub> Resolution	70 to 720 ml/min O <sub>2</sub> 1 ml/min O <sub>2</sub>	20° to 60°C 375 to 780 mm/Hg 10% to 95% non-condensing	<b>Storage Environment</b>
<b>Physical Dimensions</b>	Size (L x W x H)	20x30x10 cm	20 to 60 °C 375 to 795 mm Hg 10 to 95% RH non-condensing	<b>Part # / Product Description</b>
				9FG0080 - REEVUE Printer Kit 9FG0023 - REEVUE Basic Package 9FG0024 - Box of 20 MetaBreathers 9FG0025 - Box of 100 MetaBreathers



\*RMR is calculated using the Weir equation within assumed RQ=0.83. RMR=6.925 x VO<sub>2</sub>. Weir, J.B. New Methods for Calculating Metabolic Rate with Special Reference to Protein Metabolism. J Physiology, 1949. 109:pages 1-9

Deltatrac is a trademark of Datex-Ohmeda Corporation, Helsinki Finland

Call OPS Medical  
1(800) 430-0048

An Authorized Korn Distributor