

**MODEL 6550  
Fold-up Portable  
Wheelchair  
Scale**



**Lightweight digital  
fold-up scale for  
heavyweight tasks.**

# Portable Fold-up Scale For Special Care

This sturdy but lightweight scale offers a wide range of versatility. It can be used with wheelchairs, straight-legged chairs or as a stand-on scale, and its sloped ramps offer versatility by allowing loading from either side.

The 6550's vertical fold-up feature equals minimal floor space and allows easy transportation in four directions. Simply fold the column down, latch securely into position and it's ready to lift and move or store in an upright position. The heavy-duty wheels contribute to its ease of mobility and lack of needed storage space.

The large, easily accessed weighing platform and low angle ramps make the 6550 ideal for patient comfort and security.



The easy-to-use 758C weight display is capable of simply displaying weight or printing and sending the patient's ID, weight, BMI, time, and date to a PC. It calculates BMI (Body Mass Index) by simply entering the patient's height while weighing. The 758C also provides a crisp LCD display with battery power which adds to its flexibility.



## Specifications

Capacity:	800 lb x 0.2 lb or 300 kg x 0.1 kg (switchable)	Platform (flat):	30"W x 31.7"L (76cmW x 81cmL) with nonslip removable mat, entrance and exit ramps
Color:	White	Keypad:	19 key, color-coded, sealed membrane
Display:	Liquid crystal, reflective type	Sensitivity:	1.2V/d
Weight:	5-digit, 7-segment 0.5" high	Resolutions:	5,000 divisions, maximum
Bar Graph:	8-segment 0.3" high	Power:	6 "C" size Alkaline or Nickel-Cadium cells (not included) or optional AC adapter
BMI:	3-digit, 7-segment 0.5" high. Range: 15 - 60	Temperature Range:	+14 to +104° F/-10 to +40° C
Height ID:	11-digit entered, 4-digits displayed, 7-segment 0.3" high. Range: 3 to 7/90 to 215cm	Enclosure:	Painted steel
Dimensions (closed):	42"H x 41.5"W x 12.5"D (107cmH x 105cmW x 32cmD)	Net Weight:	61 lb (28 kg)
		Ship Weight:	91 lb (41 kg)

Cardinal Scale reserves the right to enhance, improve, or modify feature and specifications without prior notice. All registered trademarks are the property of their respective owners.



**Sold By:**

**OPS Medical, LLC** An Authorized Detecto Distributor  
**Free: 1(800) 430-0048 Local: 1(410) 553-9293**  
**email sales@opsmedical.com**  
**www.opsmedical.com**