

Cervico 2000™

Traction on the Move

Customizable Dynamic Ambulatory Traction

Cervico 2000™ is an ambulatory traction device for the treatment of acute and chronic neck pain. Easily tailored to the specific needs of each patient, Cervico 2000 is one of the few devices that enables both symmetric and asymmetric traction to the cervical spine.

Cervico 2000™ goes beyond pain relief to targeting the actual cause of the pain by actively promoting the healing of such disorders and injuries as herniated discs, torticollis and whiplash.

During treatment, patients are completely mobile, encouraged to perform light activities and exercise for active rehabilitation.

Modes of Traction:



Symmetric Traction



Asymmetric Traction

Benefits & Advantages

- Safe and effective treatment
- No drugs, hospitalization or surgery needed
- Complete freedom of movement
- Only 15 minutes per treatment session
- Compact, lightweight and user-friendly
- Reimbursable for clinic use and patient purchase
- Portable & Affordable



Compact. Customizable. Cost-Effective.

Results That Count

Clinical Facts

Traction decreases intra-discal pressure on the nerve root and dural sac. Disc disorders are accompanied by muscle spasm and reduced transfer of fluid. Due to the disc's unique structure, stationary body positioning causes reduction in fluid exchange and accelerates disc degeneration. Research shows that transfer of fluid is essential for the normal biomechanics of the disc and is possible when the spine is in motion.

The Traction on the Move concept is based on combining motion with traction which simultaneously decreases pressure, enables fluid exchange, improves Range of Motion (ROM) and promotes recovery. Cervico 2000 enables application of distraction force directly to the cervical spine alongside the cervical lordosis while retaining patient mobility.

X-Ray Results: Cervico 2000™ widens intervertebral space and releases compression



Before Treatment

The cervical spine before treatment shows hyperlordosis with nerve compression in the intervertebral foramen.



During Treatment

The cervical spine during Cervico 2000 treatment exhibits hyperflexion. The intervertebral foramen is open, there is no nerve compression and the intervertebral space is widened.

About Meditrac

Meditrac is a leader in the development, testing and manufacture of innovative medical devices for the ambulatory non-surgical treatment of spinal disorders. The "Traction on the Move" concept, developed by a team of orthopedic surgeons, combines the benefits of traction therapy with freedom of movement for rapid relief and a sense of well-being.



Call: 1(800) 430-0048

Distributed by: OPS Medical, LLC

www.opsmedical.com

Indications for Use

- Neck Pain of Discogenic Origin
- Neck Pain of Osteoarthritic or Degenerative Origin
- Whiplash Injury
- Torticollis
- Sports Injury
- Job-Related Neck Pain

Cervico 2000™ includes:

- An adjustable one-size-fits-all device with a quick release mechanism
- A carrying case, user guide, and instructional CD-ROM
- A one year warranty



Cervico 2000 carries the CE mark. Cervico 2000 is registered with the United States Food and Drug Administration. United States Federal Law permits the sale of this device (within US) on a doctor's prescription only.

