

**Material Safety Data Sheet**  
 May be used to comply with  
 OSHA's Hazard Communication Standard,  
 29 CFR 1910.1200. Standard must be  
 consulted for specific requirements.

**U.S. Department of Labor**  
**Occupational Safety and Health Administration**  
 (Non-Mandatory Form)  
 Form Approved  
 OMB No. 1218-0072



**IDENTITY** (As Used on Label and List)

SCAN® ULTRASOUND GEL

*Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.*

**Section I**

**Manufacturer's Name**

PARKER LABORATORIES, INC.

**Address (Number, Street, City, State and ZIP Code)**

286 ELDRIDGE ROAD

FAIRFIELD, NJ 07004

**Emergency Telephone Number**

**Telephone Number for Information**

(973) 276-9500

**Date Prepared**

SEPTEMBER 2006

**Signature of Preparer (optional)**

**Section II – Hazardous Ingredients/Identity Information**

| Hazardous Components (Specific Chemical Identity; Common Name(s)) | OSHA PEL | ACGIH TLV | Other Limits Recommended | % (optional) |
|---|----------|-----------|--------------------------|--------------|
|---|----------|-----------|--------------------------|--------------|

NONE

**Section III – Physical/Chemical Characteristics**

|                                |     |  |      |
|--------------------------------|-----|--|------|
| <b>Boiling Point</b>           | N/D | <b>Specific Gravity (H<sup>2</sup>O = 1)</b> | 1.01 |
| <b>Vapor Pressure (mm Hg.)</b> | N/D | <b>Melting Point</b>                         | N/D  |
| <b>Vapor Density (AIR = 1)</b> | N/D | <b>Evaporation Rate (Butyl Acetate = 1)</b>  | N/D  |

**Solubility in Water**

SOLUBLE

**Appearance and Odor**

VISCOUS, CLEAR, AQUEOUS GEL; LIGHT GREEN COLOR; PERFUMED ODOR

**Section IV – Fire and Explosion Hazard Data**

**Flash Point (Method Used)**

NON-FLAMMABLE

**Flammable Limits**

N/A

**LEL**

N/A

**UEL**

N/A

**Extinguishing Media**

USE EXTINGUISHING MEDIA APPROPRIATE FOR SURROUNDING FIRE.

**Special Fire Fighting Procedures**

NONE

**Unusual Fire and Explosion Hazards**

NONE

