

Special Offer \$6,985.99% 24 M

Eval Express™

The Chiropractor's Essential Survival Tool Kit



- ⇒ Completely portable
- ⇒ Only 12" x 6" x 6"
- ⇒ Weighs less than 3 lbs

Eval Express empowers you to

- | | |
|-----------------|--|
| Wireless | • Provide objective justification for treatment and rehabilitation |
| Range of Motion | • Show functional effects of subluxation |
| Muscle Testing | • Quantify pain sensitivity |
| Algomtry | • Substantiate efficacy of treatment |
| Objective | • Motivate continuation of care |
| Functional | • Defend yourself against post-payment audits |
| Testing | |

Over for more information ►►

Call Today for a FREE
Online Demonstration.



Phone: 800.430.0048, 410.553.0903
Fax: 410.510.1749
sales@opsmedical.com,
www.opsmedical.com



Objective Testing Made Easy

- Intuitive Software Takes Only One Hour to Learn
- Completely Wireless System with Reliable Medical Grade Technology
- Multimedia Help Videos Demonstrate Every Standard Test
- Innovative "Repeat" Button Expedites Follow-up Exams
- Free Training, Support and Software Upgrades for Six Months
- Compact Carrying Case and Complete Accessory Kit



Faster, Easier Reports

Creating reports is a snap with Eval Express' Quick-Click selection. You simply click the elements of the report you want or streamline the process with easy-to-make templates.

Report Features

- System reports include preformatted charts
- Narrative reports import into Microsoft Word
- Colorful charts and graphs
- Real time Muscle testing force curves
- Normative comparisons
- Glossary of terms
- Reusable templates

Progress Reports

Concisely demonstrate how the patient has progressed over time using text and bar graphs.

Eval Express

Range of Motion

Objectively quantify deficits in the spine and extremities and establish baselines.

- Dual and single inclinometry
- Compliant with AMA protocols
- Automatic Normative comparisons
- Test to AMA spine ROM validity criteria or to CV

Muscle Testing

Objectively quantify muscle strength, identify weakness patterns, establish baselines.

- Real-time force curves
- Select tests by region, nerve root or peripheral nerve
- Automatic bilateral comparisons and CV calculations
- Document muscle performance for all grades
- Set test thresholds and rest times

Algometry

Objectively evaluate and document pain sensitivity and tolerance.

- Built in Fibromyalgia test protocols
- Pacer for proper force application
- Automatic Normative and bilateral comparisons
- Algometry adapter with replaceable soft rubber tips



Call Today for a FREE Online Demonstration.

© JTECH Medical. October 2007.



Phone: 800.430.0048, 410.553.0903
Fax: 410.510.1749
sales@opsmedical.com,
www.opsmedical.com