

Measurement Report

12/30/01

Patient Information

First Name.....**Joan**
 Family Name.....**Smith**
 Age.....55y1m
 Height.....5'5"
 Weight.....138.0 Lb
 BMI.....23

ID.....122-12-1234
 Gender.....Female
 Menopausal age.....52
 Referring Physician.....Dr. Robinson
 Reporting Physician.....Dr. Costner
 Operator.....Lynn Moore


Measurement Results

Site..... **DISTAL RADIUS**

SOS.....4081[m/sec]
 Speed of Sound

T-Score.....-0.7
 Number of standard deviations from the mean of healthy young adults

Z-Score.....0.0
 Number of standard deviations from the mean of age matched population



Findings and Comments

- Bone strength is within the normal range.
- Bone strength is low; Medical attention is recommended.
- Bone strength shows deterioration; Medical attention is recommended.

Signature _____

Device Details

System S/N.....20092194 Probe S/N.....CMY2330
 S/W.....2.3.2645 Print Date.....12/30/01
 RDB.....American, Female

Sunlight Medical Inc

Somerset NJ, 08873, USA

Tel: 732-560-8770 Fax: 732-560-8771
 E-mail:Joanne@sunlightnet.com